

HG Series



Top tray design for accessories storage



Advanced cooling technique lowers inner temperature



Sustainable purity and flow rate for the device lifespan



Specifications

Flow Rate	0-5L/min	0-8L/min	0-10L/min
Net Weight	27/30kg	30kg	31kg
Purity	93±3%	93±3%	93±3%
Working Voltage	230V	230V	230V
Power Consumption	350W/480W	480W	680W
Outlet Pressure	0.04-0.08Mpa	0.04-0.08Mpa	0.04-0.08Mpa
Noise Level	45dB/≤50dB	≤50dB	≤55dB
Dimension	365*400*650mm	365*400*650mm	365*400*650mm
Color	Silver (Accept customization)		
OEM&ODM	Acceptable		
Standard Function	Over heat alarm.Power failure alarm. Timing function, Working hours display		
Optional Function	Low purity alarm,Nebulizer function		