

2.0ATA Lying Type Hyperbaric Chamber



Product Features

- Three-in-one control unit that combines oxygen concentrator, air compressor, and dehumidifier together has one-click operation function
- Plasma disinfection (high efficiency and low consumption)
- Provide internal/external large LCD screen, allow you to operate the equipment by yourself. Easy to control the device inside the chamber
- Multiple safety control
- Timing function (make sure the security of oxygen session)
- Automatic decompression (balance pressure in chamber)
- Manual decompression for acceleration (for emergency)
- Regulate oxygen flow (regulate oxygen flow at any time in your oxygen session)
- LCD Touchscreen, easy to operate

Specification -- Control unit

Type	Oxygen Flow	Air Flow	Size(cm)	Weight	Voltage	Power Consumption	Display screen
3-in-1 control (oxygen concentrator, compressor dehumidifier)	10L/min	120L/min	L57*W46*H70	80kg	100~240V	1300W	LCD touchable

Specification -- Chamber

Type	Pressure	Weight	Size(CM)	Material	Color	Zipper	Pressure Gauge	Frame
Soft chamber lying type	2.0ATA(14.5PSI)	14.8kg	D80*L200	TPU	White	2	inside and outside	Internal metalframe